Instruction Of Use

AMERWAY OXYGEN REGULATOR: MODEL: IVR-25L-O (ATYF – 10)

Device description:

Intended Use:

The device that adjust the output gas pressure, flow rate, oxygen concentration and other parameters according to a predetermined. For oxygen emergency use and for oxygen-absent patient breathing use.

Indications

1. Respiratory diseases: pneumonia, bronchitis, chronic bronchitis, viral respiratory infections, asthma, emphysema, cor pulmonale, etc.
2. Cardiovascular and cerebrovascular diseases: hypertension, heart disease, coronary heart disease, myocardial infarction, cerebral thrombosis, cerebral ischemia, cerebral vertigo, atherosclerosis, etc.
3. Plateau hypoxia disease: plateau pulmonary edema, acute mountain sickness, chronic mountain sickness, plateau coma, plateau hypoxia, etc.
4. People vulnerable to hypoxia: middle-aged and old people, pregnant women, students who have been engaged in mental work for a long time, company staff and office cadres; People who have been engaged in underground work or confined space for a long time, and those who have excessive exercise, fatigue and insufficient oxygen supply, etc.
5. Other people in need of oxygen therapy: weak and sickly people with poor immunity, heat stroke, gas poisoning, drug poisoning, etc.

Oxygen Regulator Components:

1. Flow knob (turn on or off outlet and adjustable the flow rate): FLOW RATE 0-25L
2. Gas Mouth: Filling connector (filling oxygen in the cylinder): FILLING CONNECTOR: CGA540
3. The outlet (quick socket): QUICK CONNECTOR: DISS
4. Connected with the cylinder thread: CYLINDER CONNECTOR: PZ27.8
5. Pressure gauge: PRESSURE GAUGE: 0-3000 Psi
6. Main switch knob: ON/OFF KNOB

Installation

- (1) Installed oxygen regulator installed before install bottle, please make sure the seal being loaded on the thread, sealing ring or teflon wrap, please check whether the damaged, if damaged, please replaced in a timely manner.
- (2) Take off the regulator handle and use a suitable tool for valve assembling. Valving square is located above inlet connection. Fit the handle with the screw.
- (3) To assembly, please follow the tightening torques given by EN ISO 13341 standards or recommended by the supplier of cylinders.

Operating instructions:

- (1) General use: Open main switch knob, Fit the mask or nasal, turn on the flow knob, adjustable the flow rate.
- (2) High traffic usage: Open main switch knob, use a hose with suitable connector fit the outlet (quick socket), it is about 100psi flow rate oxygen.
- (3) Filling the gas: Open main switch knob, then use a suitable hose connect the Filling connect to fill the cylinder.
- (4) Use the oxygen mask and connecting pipe connection is good, and the oxygen tube at one end connected to the air outlet end on portable valve, rotating open black knob to the right, the pressure gauge pointer rotation, according to the cylinder pressure (chart), valve showed normal work, had oxygen output and oxygen mask to oxygen mask wear well, turn the knob, choose needed flow. To turn off, please turn black knob to the left to the end, and to ensure that no oxygen both senders.
- (5) Inflation: Directly connect the inflation pipe to the inflation port of the portable valve, and connect the other end to the inflation station equipment or other inflation equipment. Note that when loosening the oxygen bridge from the inflation port, it is best to use a wrench to prevent the inflation port from loosening and causing air leakage.

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