

## ANAESTHETIC GAS SCAVENGING SYSTEM






### Features




- GENTEC® Anaesthetic Gas Scavenging System is used for the safe removal of excess anaesthetic gas and anaesthetic vapors
- Color-coded for easy gas service identification according to ISO 32
- Complies with ISO 7396-2 and ISO 9170-2
- Flow: 50 l/min, 80 l/min (Pressure drop: 0.15~0.29 psi)
- Air driving flow : 10~15 l/min (Pressure: 58 psi)
- 100% tested and individually packaged

## Dimensions (mm)




3816E-AGSS-B



3826E-AGSS-F



3826E-AGSS-W



## Ordering Information

| Part Number | Model Number          | Maximum Flow | Exhaust Diameter | Installation  |
|-------------|-----------------------|--------------|------------------|---------------|
| A65216090   | 3816E-AGSS(1L)-HB16-B | 50 L/min     | I.D 16 mm Hose   | Console Mount |
| A65216100   | 3816E-AGSS(1H)-HB16-B | 80 L/min     | I.D 16 mm Hose   | Console Mount |
| A65226010   | 3826E-AGSS(1L)-HB22-W | 50 L/min     | I.D 22 mm Hose   | Box Mount     |
| A65226020   | 3826E-AGSS(1H)-HB22-W | 80 L/min     | I.D 22 mm Hose   | Box Mount     |
| A65226030   | 3826E-AGSS(1L)-HB16-F | 50 L/min     | I.D 16 mm Hose   | Panel Mount   |
| A65226040   | 3826E-AGSS(1H)-HB16-F | 80 L/min     | I.D 16 mm Hose   | Panel Mount   |