

SUCTION THERAPY

Suction Regulator

Solutions for Life



880VS Series Venturi Suction Regulator

Glow in the Dark Gauge

- 2" gauge, easy to read
- Color coded flow ranges
- Dual Scale mmHg and kPa

Safety Device

- Positive pressure relief valve provides safe use
- ISO 10079-3 Compliant

Control Knob

- Easy adjustment of vacuum pressure
- No Leakage

Suction Trap

- Prevents contamination caused by overflow

Type of Adapter

- None: Without adapter
- OH: Ohmeda® adapter
- DS: DISS hex adapter
- DH: DISS Handwheel adapter
- CH: Chemetron® adapter
- PB: Puritan-Bennett® adapter
- FS: French (NF S 90-116) adapter
- GS: German (DIN 13260-2) adapter
- BS: British (BS5682-1998) adapter
- JIS: Japanese Style adapter
- SIS: Australian (AS2896) adapter

Venturi Suction

- Delivers steady and high vacuum flow from a medical air or oxygen gas source
- Control Range: -1 kPa ~ -80 kPa

Suction Connector

- GENTEC® provides a wide range of fittings to meet customers' needs



880VS-760-00-A-T

Technical Specification

Dimension (H x W x D)	90 x 90 x 60 mm
Weight	490 g
Range	-1 kPa ~ -80 kPa
Maximum Suction	> 25 L / min
Source of Pressure	300 ~ 500 kPa
Vacuum Control	Needle valve control venturi
Vacuum Gauge Range	0 ~ -100 kPa / 0 ~ 760 mmHg

Ordering Information

Ex: 880VS	- 760	- 00	- A	- T
Series	Measuring Range	Vacuum Connection	Driving Gas	Suction Connection
880VS	760: 0~760 mmHg	00: No adapter, 1/8" NPT (F) OH: 3111U (Ohmeda) adapter DH: 3222U (DISS hand-tight) adapter DS: 3102 (DISS Hex) adapter CH: 3103 (Chemetron) adapter PB: 3104 (Puritan-Bennett) adapter FS: 3105 (FS) adapter GS: 3116 (GS) adapter BS: 3117 (BS) adapter JIS: 3108 (JIS) adapter SIS: 3109 (SIS) adapter	A: Medical Air O: Oxygen	00: No adapter, 1/8" NPT (F) T: With suction trap assembly (complete) TD: With suction trap (complete DISS Female nut) TH: With suction trap (complete DISS female vacuum hand-tight) H: With tubing nipple (5/16" hose barb) D: With DISS male vacuum N: With NIST adapter F: With trap fitting only (no trap)